



Introducing the new

# ENFit connection system

The new standard in enteral feeding connections

## The ABC of ENFit Transition Connectors


Ensuring a smooth transition to the ENFit Connection System



All three transition adapters (A, B, and C) will be supplied in a **single pack** included within the Feeding Set packaging. Adapters will only be needed during the transition period. From March 2016, new ENFit syringes and feeding tubes will be available to connect directly to the ENFit giving set. Adapters will be available until December 2016.

### Post-transition connection

ENFit connection only (no adapters required) once all feeding tubes and syringes are transitioned to the ENFit System



For more information, call the **Hospital to Home Helpline** on **Freephone 1800 22 11 66** or visit **www.abbottnutrition.ie**